

Mölnlycke® dressing cutting guide



Tip of toe
10x10 cm







Hand
10x10 cm
or 10x20 cm



Around tubes
10x10 cm
or 10x20 cm



Elbow
10x10 cm



Ear
10x10 cm






Between toes
10x10 cm



Heel
10x10 cm



General all purpose
for younger patients
10x10 cm



Proving it every day

At Mölnlycke®, we deliver innovative solutions for managing wounds, improving surgical safety and efficiency and preventing pressure ulcers. Solutions that help achieve better outcomes and are backed by clinical and health-economics evidence.

In everything we do, we are guided by a single purpose: to help healthcare professionals perform at their best. And we're committed to proving it every day.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00. The Mölnlycke, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2020 Mölnlycke Health Care AB. All rights reserved. XXXXXXXX

