



Tubifast® 2-Way Stretch®: What size to choose?


RED LINE




GREEN LINE




BLUE LINE



YELLOW LINE



PURPLE LINE



Limb	Limb circumference	Product line	Product length
Small limbs	9 – 18 cm 3.5 – 7 in	RED LINE	1 m – 10 m
Small & medium limbs	14 – 24 cm 5.5 – 9.5 in	GREEN LINE	1 m – 3 m – 5 m – 10 m
Large limbs	24 – 40 cm 9.5 – 16 in	BLUE LINE	1 m – 3 m – 5 m – 10 m
Extra-large limbs, heads, children's trunks	35 – 64 cm 14 – 25 in	YELLOW LINE	1 m – 3 m – 5 m – 10 m
Adult trunks	64 – 130 cm 25 – 51 in	PURPLE LINE	1 m – 5 m – 10 m

Find Tubifast® 2-Way Stretch® application videos at <https://www.molnlycke.com/education/wound-areas/wound-healing/cutting-guide-for-dressings/>

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00. The Mölnlycke, Tubifast and TwoWay Stretch trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2020 Mölnlycke Health Care AB. All rights reserved.

